




***Surgical Table
Attachments & Accessories***

Product Catalogue



› +61 3 9764 0900
sales@taylorsurgical.com.au
www.taylorsurgical.com.au

TABLE OF CONTENTS



UPPER LIMB SUPPORT

TV Arm boards.....	4
TV 'J' Board with Pad.....	5
TV Carter Braine Arm Support.....	6
TV Elbow Arm Support.....	7
TV Adult Hand & Arm Table.....	8
TV Standard Arm board With Pad.....	9
TV Micro Surgical Arm & Hand Table.....	10

ANAESTHETIC

TV Anaesthetic Support Screen.....	11
TV Anaesthetic Tube Holder.....	12
TV Hinged "Fan Out" Anaesthetic Screen.....	13
TV Transducer Holder.....	14

LOWER LIMB SUPPORT

TV Latchlock Light Weight Board.....	15
TV Nissen Seat & Light Weight Board.....	16
TV Peri Tray Attachment.....	17
TV Popliteal Support.....	18
TV Lower Limb Prep Support.....	19

LATERAL & HIP SUPPORT

TV Curved Lateral Support.....	20
TV Lateral Hip Positioner.....	21

HEAD SUPPORT

TV Narrow Headrest.....	22
TV 4085 Neuro Head Rest Adaptor.....	23

NEONATE SUPPORT

TV Neonate Infant Table Attachment.....	24
---	----

CLAMPS

TV Multipurpose Clamp.....	25
TV Clark Socket Clamps.....	26
TV Gantry Shoulder Clamp.....	26

SIDE RAILS

TV 4085 'Denyer' Side Rails.....	27
TV 4085 'Eschmann' Side Rails.....	27
TV 3085 'Denyer' Side Rails.....	28
TV 3085 'Eschmann' Side Rails.....	28
TV Side Rail Extenders.....	29
TV Side Rail Extended Length.....	29

ACCESSORIES

TV Accessory Storage Cart.....	30
TV Hand Table Telescopic Adjustable Leg.....	30
TV 2 Piece Velcro Restraint Strap with Side Rail Mount ...	31

TV ARMBOARDS

KEY FEATURES

- The radiolucent carbon fibre armboard is easily attached or removed by a gravity latch lever from the surgical table rail
- The positioning component is fully enclosed within the machined housing protecting against denting, breakage and damage when dropped
- Includes a (762mm) Velcro arm restraint strap



PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-DY43-100	TV Steris Armboard (Denyer Rail)
DC-ES43-100	TV Steris Armboard (Eschmann Rail)



UPPER LIMB SUPPORT

TV 'J' BOARD WITH PAD

KEY FEATURES

- The Lightweight "J" curved arm support (100mm x 350mm) provides forearm, elbow or wrist support throughout Anaesthesia or surgical procedures.
- Comes with protective pad

PRODUCT ORDER DETAILS

CAT. No.	DESCRIPTION
DC-JBOARDM	TV "J" Board with Pad





TV CARTER BRAINE ARM SUPPORT

KEY FEATURES

- Â The arm support is appropriate for unilateral prone, supine and seated positioning of patients
- Â The universal ball joint with concaved forearm construction allows complete flexibility when positioning the patients arm

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV310	TV Carter Braine Arm Support

OTHER ATTACHMENTS REQUIRED (Dependant on Rail Size)

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

* Ordered separately



TV ELBOW ARM SUPPORT

KEY FEATURES

- The elbow arm support is appropriate for unilateral prone, supine and seated positioning of patients and most commonly used for "elbow arthroscopy"
- The universal ball joint with concaved forearm construction allows complete flexibility when positioning the patients arm

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV311E	TV Elbow Arm Support



OTHER ATTACHMENTS REQUIRED (Dependant on Rail Size)

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

*Ordered separately

TV ADULT HAND & ARM TABLE

KEY FEATURES

- Â The radiolucent (400mm x 800mm) adult hand & arm table provides unobstructed C-Arm and surgeon access for major and minor surgical procedures
- Â The size also accommodates minor Vascular or Orthopaedic procedures
- Â Can have "optional telescopic adjustable leg" fitted to it for extra support



PRODUCT DETAILS

CAT. No. DESCRIPTION

CAT. No.	DESCRIPTION
DC-TV43-500	TV Adult Hand & Arm Table

OTHER ATTACHMENTS REQUIRED

Dependant on Rail Size

CAT. No. DESCRIPTION

DC-US600A-DE	TV multi-purpose clamp (US Rail)
--------------	----------------------------------

DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
--------------	--------------------------------------

DC-DY605S	Clark Socket (Denyer Rail)
-----------	----------------------------

DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
--------------	--

DC-ES605S*	Clark Socket (Eschmann Rail)
------------	------------------------------

* Ordered separately



TV STANDARD ARM BOARD WITH PAD

KEY FEATURES

- The single pin radiolucent (280mm x 750mm) standard length arm board provides the platform for Anaesthesia cannulation, unobstructed C-Arm and surgeon access for minor surgical procedures
- The size also accommodates minor Vascular or Orthopaedic procedures

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-STD43-100	TV Standard Arm Board with Pad



OTHER ATTACHMENTS REQUIRED (Dependant on Rail Size)

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

* Ordered separately



UPPER LIMB SUPPORT

TV MICRO SURGICAL ARM & HAND TABLE

KEY FEATURES

- The "Hourglass" shaped (650mm x 850mm) Micro surgical arm & hand table provides stability for the surgeon when performing complex Vascular and Orthopaedic procedures
- The radiolucent construction allows C-Arm access of the extremities and comes with an adjustable telescopic leg fitted for extra support



PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV43-600	TV Adult Hand & Arm Table

OTHER ATTACHMENTS REQUIRED *Dependant on Rail Size*

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

** Ordered separately*



TV ANAESTHETIC SUPPORT SCREEN

KEY FEATURES:

- Â The anaesthetic support screen is designed to support surgical drapes over the patients head so the anaesthetist can maintain full access clearance to the face
- Â The length of the "L shaped" bar allows adjustable height flexibility



PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-ANASCREEN	Anaesthetic Support Screen

OTHER ATTACHMENTS REQUIRED *Dependant on Rail Size*

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

** Ordered separately*

ANAESTHETIC




TV ANAESTHETIC TUBE HOLDER

KEY FEATURES:

- The tube holder is purposely design to support the airway management tubing from the anaesthetic machine throughout surgery

PRODUCT ORDER DETAILS:


CAT. No.	DESCRIPTION
DC-TV TUBEH	Anaesthetic TubeHolder



TV HINGED “FAN OUT” ANAESTHETIC SCREEN

KEY FEATURES

- The hinged (fan out) anaesthetic screen attaches to the neonate infant table supporting surgical drapes over the patients head whilst allowing the Anaesthetist to raise the screen upwards to maintain full airway and patient access
- Can be positioned in three (3) 100mm increments dependant on the infant’s length



PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV31-565	Hinged “fan out” anaesthetic screen

OTHER ATTACHMENTS REQUIRED

CAT. No.	DESCRIPTION
DC-TV32-400	Neonate infant table attachment

TV TRANSDUCER HOLDER

KEY FEATURES

- The Transducer Holder provides a stable post for securing a transducer at required height for Neuro and Cardiac surgery
- The curved design makes it easy for Anaesthesiologist to position and view the transducer

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV26-300	TV Transducer Holder

OTHER ATTACHMENTS REQUIRED

Dependant on Rail Size

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

**Ordered separately*



TV LATCHLOCK LIGHT WEIGHT BOARD

KEY FEATURES

- Â The specially designed "Latchlock" light weight board attaches to the leg section insert tabs of a STERIS C-Max, 4085 / 5085 surgical table
- Â The unique design, shape & length of the latchlock board caters for patients going into stirrups allowing full support of the legs and lower limbs
- Â Most commonly used in Urology, Gynaecology, General & Colorectal procedures
- Â Comes complete with a 2 Inch support pad



PRODUCT ORDER DETAILS

CAT. No.	DESCRIPTION
DC-LATCHLOCK	Latchlock light weight board



TV NISSEN SEAT & LIGHT WEIGHT BOARD

KEY FEATURES:

- Â The specially designed "Bariatric" Nissen Seat attaches to the leg section insert tabs of a Steris C-Max / 4085 / 5085" surgical table
- Â Designed around leading Gastric Band surgeons, the unique shape contours around the patient's thighs reducing the risk of pressure point
- Â The nissen seat comes with its own light weight board for leg and lower limb support
- Â It comes complete with a 2 Inch support pad

PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-NISSENSEAT	Nissen seat and light weight board



Â Nissen Seat



Â Light Weight Board



Â Nissen Seat (Sideview)



TV PERI TRAY ATTACHMENT

PRODUCT DETAILS

KEY FEATURES

- Â The Peri Tray attaches to the leg section insert tabs of a "C-Max /4085 / 5085" surgical table
- Â The Peri Tray is designed as an attachable instrumentation tray with the option of being repositioned or articulated via the tables electronic hand control
- Â Most commonly used in Urology, Gynaecology, General & Colorectal procedures.
- Â Comes with a removable tray for easy cleaning or sterilization

CAT. No.	DESCRIPTION
DC-PERITRAY	Peri Tray Attachment



TV POPLITEAL SUPPORT

KEY FEATURES

- The popliteal support is appropriate for joint replacement, general, unilateral prone, supine and seated patient positioning

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-DY300	Popliteal Support

OTHER ATTACHMENTS REQUIRED *Dependant on Rail Size*

CAT. No.	DESCRIPTION
DC-US600A-DE	US Rail
ST-BF08-300	Clark Socket US Rail
DC-DY600A-DE	Denyer Rail
DC-DY605S	Clark Socket Denyer Rail
DC-ES600A-DE	Eschmann Rail
DC-ES605S*	Clark Socket Eschmann Rail

**Ordered separately*



LOWER LIMB SUPPORT



TV LOWER LIMB PREP SUPPORT

KEY FEATURES

- The Lower Limb support is appropriate for unilateral prone, supine and seated positioning of patients limbs for surgical skin preparation

PRODUCT DETAILS


CAT. No.	DESCRIPTION
DC-625S	TV Lower Limb Prep Support

OTHER ATTACHMENTS REQUIRED

Dependant on Rail Size

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
ST-BF08-300	Clark Socket (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

**Ordered separately*



LATERAL & HIP SUPPORT

TV CURVED LATERAL SUPPORT

KEY FEATURES:

- The TV curved lateral support allows you to position & restraint patients throughout surgical procedures
- Adjusting both vertically and horizontally it's curved shape contours to the patients – commonly used in Orthopaedic, General, Bariatric etc...

PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-ES305S	TV Curved Lateral Support - Small Size
DC-ES305M	TV Curved Lateral Support - Medium Size

OTHER ATTACHMENTS REQUIRED

Dependant on Rail Size

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)


**Ordered separately*



• TV Curved lateral support



• TV Curved lateral support showing patient



• TV Curved lateral support attached to table

TV LATERAL HIP POSITIONER

KEY FEATURES

- Â The TV lateral Hip Positioner support allows you to position & restraint patients throughout surgical procedures
- Â Adjusting both vertically and horizontally the round disc pads can be applied to the preferred location allowing surgical site access – commonly used in Orthopaedic, General, Bariatric etc...

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV315	TV Lateral Hip Positioner



OTHER ATTACHMENTS REQUIRED *Dependant on Rail Size*

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
ST-BF08-300	Clark Socket (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyer Rail)
DC-DY605S	Clark Socket (Denyer Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)
DC-ES605S*	Clark Socket (Eschmann Rail)

*Ordered separately



HEAD SUPPORT




TV NARROW HEADREST

KEY FEATURES

- The Narrow Headrest is ideal for Ophthalmic, ENT, Neurological General & Plastic procedures
- The headrest inserts to the back section of a Steris "C-Max/ 3085 / 4085 / 5085" surgical table
- Designed around the different aspects of surgery allowing procedural access for surgeons
- Comes complete with a 2" Inch support pad

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV59-200	TV Narrow Headrest



HEAD SUPPORT

TV 4085 NEURO HEAD REST ADAPTOR

KEY FEATURES

- The custom design adaptor accommodates new and older style Neuro Mayfield head rest, eye & ENT support systems

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV58-900	TV 4085 Neuro Head Rest Adaptor



OTHER ATTACHMENTS REQUIRED

Mayfield Head Rest Support - Sold Separately



TV NEONATE INFANT TABLE ATTACHMENT

KEY FEATURES

- The Neonate Infant Table (400mm wide x 800mm long) securely attaches to the end of back/seat section bringing the surgeon closer to improve surgical access for Infant patients
- Comes complete with a 2 Inch support pad

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV32-400	Neonate Infant Table Attachment



CLAMPS

TV MULTIPURPOSE CLAMP

KEY FEATURES

- The TV Denyer / Eschmann multipurpose clamp allows you to support round, oval, blade or square pin table accessories
- The contoured insert cushion allows adjustable height flexibility with mid rail access to US Rails

PRODUCT ORDER DETAILS

CAT. No.	DESCRIPTION
DC-US600A-DE	TV multi-purpose clamp (US Rail)
DC-DY600A-DE	TV multi-purpose clamp (Denyers Rail)
DC-ES600A-DE	TV multi-purpose clamp (Eschmann Rail)



TV CLARK SOCKET CLAMPS

KEY FEATURES

- The TV Denyer / Eschmann Clark Socket clamp allows you to support round or oval pin accessories anywhere along the table side rail
- The sockets provide angular control and positive locking for any accessory in a plane parallel to the side rail and allows adjustable height flexibility with mid rail access to US Rails

PRODUCT DETAILS

CAT. No.	DESCRIPTION
ST-BF08-300	Steris Clark Clamp (US Rail)
DC-DY605S	TV Denyer Clark Clamp (Denyer Rail)
DC-ES605S	TV Eschmann Clark Clamp (Eschmann Rail)



TV GANTRY SHOULDER CLAMP

KEY FEATURES

- The TV Gantry clamp allows you to support round or oval pin shoulder arthroscopy traction systems
- The contoured insert cushion allows adjustable height flexibility with mid rail access to US Rails

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-US600A-GANTRY	US Rail

OTHER ATTACHMENTS REQUIRED

CAT. No.	DESCRIPTION
DC-DY600A-GANTRY	Denyer Rail
ES600A-GANTRY	Eschmann Rail



SIDE RAILS



TV 4085 'DENYER' SIDE RAILS

KEY FEATURES

- The Denyer "4085 side rail" has been designed to perfectly match the Steris 4085 Surgical Table and can be ordered with or without the popular "mid mount" rail access similar to the Steris US side rail
- The 316 high polished (6mm x 38mm) stainless steel side rails come in a set /8

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-4085-Denyer	4085 "Denyer" Side Rail Suited Steris 4085 Surgical Table
DC-4085-Denyer-MID	4085 "Denyer" Side Rail Suited Steris 4085 Surgical Table

TV 4085 'ESCHMANN' SIDE RAILS

KEY FEATURES

- The Eschmann "4085 side rail" has been designed to perfectly match the Steris 4085 Surgical Table and can be ordered with or without the popular "mid mount" rail access similar to the Steris side rail
- The 316 high polished (6mm x 31mm) stainless steel side rails come in a set /8

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-4085-Eschmann	4085 "Eschmann" Side Rail Suited Steris 4085 Surgical Table
DC-4085-Eschmann-MID	4085 "Eschmann" Side Rail Suited Steris 4085 Surgical Table



SIDE RAILS



TV 3085 'DENYER' SIDE RAILS

KEY FEATURES

- Â The Denyer "3085 side rail" has been designed to perfectly match the Steris 3085 Surgical Table and can be ordered with or without the popular "mid mount" rail access similar to the Steris US side rail
- Â The 316 high polished (6mm x 38mm) stainless steel side rails come in a set /8

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-3085-Denyer	3085 "Denyer" Side Rail Suited Steris 3085 Surgical Table
DC-3085-Denyer-MID	3085 "Denyer" Side Rail Suited Steris 3085 Surgical Table

TV 3085 'ESCHMANN' SIDE RAILS

KEY FEATURES

- Â The Eschmann "3085 side rail" has been designed to perfectly match the Steris 3085 Surgical Table and can be ordered with or without the popular "mid mount" rail access similar to the Steris side rail
- Â The 316 high polished (6mm x 31mm) stainless steel side rails come in a set /8
- Â Optional custom design

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-3085-Eschmann	3085 "Eschmann" Side Rail Suited Steris 3085 Surgical Table
DC-3085-Eschmann-MID	3085 "Eschmann" Side Rail Suited Steris 3085 Surgical Table



SIDE RAILS


TV SIDE RAIL EXTENDERS

KEY FEATURES:

- Â The TV Heavy Duty Side Rail Extender allows you to extend the side rail past the patients head
- Â Constructed from high polished 316 stainless steel, the extender rail is most commonly used for arm board attachment, side rail clamps etc...
- Â The rail is provided at a standard length of 275mm long whilst the thickness and height will be determined by the type of surgical table side rail selected

PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-TV45-600	TV US H/D Side Rail Extender (USRail)
DC-TV45-600-DY	TV Denyer H/D Side Rail Extender (Denyer Rail)
DC-TV45-600-ES	TV Eschmann H/D Side Rail Extender (Eschmann Rail)




TV SIDE RAIL (HEAVY DUTY) EXTENDED LENGTH

KEY FEATURES

- Â The TV Heavy Duty Side Rail Extender allows you to extend the side rail past the patients head
- Â Constructed from high polished 316 stainless steel, the extender rail is most commonly used for arm board attachment, side rail clamps etc...
- Â The rail is provided at an extended length of 500mm long whilst the thickness and height will be determined by the type of surgical table side rail selected

PRODUCT DETAILS

CAT. No.	DESCRIPTION
DC-TV45-600-USL	TV US H/D Side Rail Extended Length (USRail)
DC-TV45-600-DYL	TV Denyer H/D Side Rail Extended Length (Denyer Rail)
DC-TV45-600-ESL	TV Eschmann H/D Side Rail Extended Length (Eschmann Rail)



STORAGE

TV ACCESSORY STORAGE CART

KEY FEATURES:

- Â The "Unique" design provides durability with the flexibility to store your table and associated accessories safely with a touch of class. Designed to store 8 x side rail clamps, 2 x Curved Lateral Supports, 1 x Lateral Hip Positioner, 2 x Arm Boards, Latchlock or light weight board & pair of powerlift stirrups
- Â The cart can be ordered to suit US, Denyer, Maquet, Eschmann or custom fit side rail accessories

PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-TV46-308	TV Accessory Storage Cart


TV HAND TABLE TELESCOPIC ADJUSTABLE LEG

KEY FEATURES:

- Â The telescopic adjustable leg was designed for the "Hourglass" shaped (650mm x 850mm) Micro surgical Arm & Hand Table to provide stability for the surgeon when performing complex Vascular and Orthopaedic procedures
- Â It can be fitted to the adult & hand table or retro fitted to any existing hand table the health facility made have

PRODUCT DETAILS:

CAT. No.	DESCRIPTION
DC-TV67-308	TV Hand Table Telescopic Adjustable Leg



Â Currently under design

TV 2 PIECE VELCRO RESTRAINT STRAP WITH SIDE RAIL MOUNT

KEY FEATURES:

- The 2 piece velcro patient restraint strap has been designed to perfectly match any of the (US Rail) and all of the Steris Surgical Table fleet
- Comes standard with their popular "Mid Mount" rail access

PRODUCT ORDER DETAILS:

CAT. No.	DESCRIPTION
DC-US620	TV 2 Piece Velcro Restraint Strap with Side Rail Mount

(Eschmann / Denyer / Maquet – Option also available)



Proudly
manufactured by

